

Philips Medical Systems

P.O. Box 10000, 5680 DA Best, The Netherlands

XJB 148 – 7107
Tel: +31402763274
Fax: +31402763085

TO WHOM IT MAY CONCERN

November 13, 2005.

Iradimed Corp., Winter Park, Florida, USA, market the mRidium 3850 MR Compatible IV pump as compatible with the Philips MRI-systems.

Philips Medical Systems Nederland B.V., at their facilities in Best, The Netherlands has performed limited testing on one sample mRidium 3850 in October 18, 2005.

The mRidium 3850 IV pump was unit: model 3850, sn0025

The battery Charger power supply was unit: model Charger 1120, sn0028

Philips found that the IV pump did not interfere with the operation of Philips MRI systems Achieva 1.5T, Achieva 3.0T Nova and the Panorama 1.0T. When the IV pump was operated through the battery charge power supply this did not result in any interference with the Philips MRI systems Achieva 1.5T, Achieva 3.0T Nova. Considering the similarity between the tested MRI-systems and other Philips MRI-configurations, Philips expects that the injector operated through the battery charge power supply will not interfere with the operation of the Philips MRI-systems Gyroscan NT 1.5T, Intera 1.5T and 3.0T the Achieva 1.5T and 3.0T.

The results of this testing may not apply to other samples of this model IV pump nor to the IV pump operated through the AC mains power supply, and you may wish to contact Iradimed to verify the applicability of this testing to the particular IV pump that you are considering.

Ref: VR test report, PMS, XJR-153-4686
test report Iradimed, received on, October 28, 2005.

Dr. Hans Engels

MR Safety Director
Philips Medical Systems Nederland B.V
The Netherlands

