

MEDICATIONS	PRIMARY INFUSIONS										BOLUS INFUSIONS										
	Dose Values						CONC Values				Dose Values						Time Values				
	Max Hard Limit	High Soft Limit	DOSE Initial Value	Dose Units	Lower Soft Limit	Min Hard Limit	Drug CONC	CONC Units	Diluent Volume	mL(s)	Max Hard Limit	High Soft Limit	BOLUS Initial Value	Dose Units	Lower Soft Limit	Min Hard Limit	Max Hard Limit	High Soft Limit	BOLUS Delivery Time	Lower Soft Limit	Min Hard Limit
amiodarone	1.5	1		mg/min	0.5		360	mg	200	mL(s)											
argatroban	10.1	10		mcg/kg/min	0.1		250	mg	250	mL(s)											
bivalirudin250mg/500	2	1.76		mg/kg/hr	0.1		250	mg	500	mL(s)											
bivalirudin 250mg/50	2	1.76		mg/kg/hr	0.1		250	mg	50	mL(s)											
bumetanide		2		mg/hr	0.1		12.5	mg	250	mL(s)	4		mg	0.25				2			1
cisatracurium	11	8		mcg/kg/min	0.5		200	mg	100	mL(s)											
dexmedetomidine	2	1		mcg/kg/hr	0.05		400	mcg	100	mL(s)											
diltiazEM	25	15		mg/hr	0.5		100	mg	100	mL(s)											
DOBUTamine		20		mcg/kg/min	1		250	mg	250	mL(s)											
DOBUTamine DS		20		mcg/kg/min	2		500	mg	250	mL(s)											
DOPamine	50	20		mcg/kg/min	1		800	mg	500	mL(s)											
DOPamine DS	50	20		mcg/kg/min	1		800	mg	250	mL(s)											
EPINEPHrine		10		mcg/min	1		5	mg	250	mL(s)											
EPINEPHrine DS		10		mcg/min	1		10	mg	250	mL(s)											
eptifibatide	2.4	2.2		mcg/kg/min	0.5		75	mg	100	mL(s)											
esmolol		300		mcg/kg/min	25		2500	mg	250	mL(s)											
fentaNYL		200		mcg/hr	0.1		1000	mcg	100	mL(s)	100		mcg	12.5				2			1
furosemide		40		mg/hr	1		250	mg	250	mL(s)											
heparin	4500	2200		Units/hr	200		25000	Units	250	mL(s)											
HYDRomorphone		4		mg/hr	0.1		50	mg	50	mL(s)	2		mg	0.01				3			2
insulin REGULAR	40	20		Units/hr	1		100	Units	100	mL(s)											
isoproterenol		10		mcg/min	1		2000	mcg	500	mL(s)											
ketamine infusion	5			mg/kg/hr	0.05		1000	mg	100	mL(s)											
labetalol		2.5		mg/min	0.05		250	mg	250	mL(s)											
lidocaine mcg/kg/min		21		mcg/kg/min	9		2000	mg	500	mL(s)											
lidocaine mg/min		5		mg/min	0.5		2000	mg	500	mL(s)											
LORazepam	19	10		mg/hr	0.25		50	mg	100	mL(s)	10		mg	1				2			
midazolam	10	7		mg/hr	0.25		25	mg	250	mL(s)	5		mg	1				2			
midazolam DS	10	7		mg/hr	0.25		50	mg	250	mL(s)	5		mg	1				2			
milrinone		0.75		mcg/kg/min	0.13		20	mg	100	mL(s)											
morphINE		15		mg/hr	0.1		50	mg	50	mL(s)	10		mg	0.25				4			
naloxone		2000		mcg/hr	250		1000	mcg	250	mL(s)											
niCARDipine		15		mg/hr	0.5		25	mg	250	mL(s)											
niCARDipine DS		15		mg/hr	0.5		50	mg	250	mL(s)											
nitroglycerin		200		mcg/min	1		100	mg	250	mL(s)											
nitroprusside		10		mcg/kg/min	0.1		50	mg	100	mL(s)											
norepinephrine		210		mcg/min	0.5		8	mg	500	mL(s)											
norepinephrine DS		210		mcg/min	0.5		16	mg	500	mL(s)											
octreotide		75		mcg/hr	25		500	mcg	100	mL(s)											
pantoprazole		9		mg/hr	7		80	mg	250	mL(s)											
phenylephrine		360		mcg/min	1		20	mg	250	mL(s)											
phenylephrine DS		360		mcg/min	1		40	mg	250	mL(s)											
procainamide		4		mg/min	1		2000	mg	500	mL(s)											
propofol	130	100		mcg/kg/min	5		1000	mg	100	mL(s)	20		mg	10				1			
rocuronium 100mg/200		15		mcg/kg/min	4		100	mg	200	mL(s)											
rocuronium 500mg/100		15		mcg/kg/min	4		500	mg	100	mL(s)											
vasopressin		6		Units/hr			40	Units	200	mL(s)											
vecuronium	2.5	1.6		mcg/kg/min	0.01		20	mg	100	mL(s)											
ANES Propofol	250			mcg/kg/min			1000	mg	100	mL(s)	50		mg					0.5			